# Tracheotomie Symposium

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
# Welkom

Lieselotte Schiltz | Business Unit Manager Atos Medical  
Benelux

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
11|20 TRACHEOTOMIE  
MEI|23 SYMPOSIUM


## WHAT IS SOCIAL CONNECTION?

The subjective experience of feeling close to and a sense of belongingness with others.

## THE BENEFITS OF HIGH SOCIAL CONNECTION:

50% increased chance of longevity

stronger gene expression for immunity (research by Steve Cole, UCLA)

lower rates of anxiety and depression

higher self-esteem and empathy

better emotion regulation skills

Social connection creates a positive feedback loop of social, emotional, and physical well being.

## THE DANGERS OF LOW SOCIAL CONNECTION

worse for health than smoking, high blood pressure or obesity

higher inflammation at the cellular level

higher susceptibility to anxiety and depression

slower recovery from disease

increased antisocial behavior and violence

suicide

Connectedness & Health: The Science of Social Connection, Emma Sappale

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







# Indicatie

# Timing tracheostomie

Predicting the need for  
admitted to an intensive care unit: A

Cristian Papuzinski<sup>a</sup>, Miguel Durante<sup>a</sup>, Catalina Tobar, MD<sup>b</sup>,  
Felipe Martinez, MD<sup>c,\*</sup>, Eduardo Labarca, MD<sup>d</sup>

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ELSEVIER

Review

Systematic review of international guidelines for tracheostomy in COVID-19 patients

Carlos M. Ch...  
Nicola...

Rome R. Lechien<sup>a,b,1</sup>, Christian Calvo-F...  
Shazia Peer<sup>a,8</sup>, Jon A. Sistiaga-Sua...  
Cammaroto<sup>a,h</sup>, Miguel M...

Tracheostomy  
in Patients  
Ventilation

What Is the Optimal Timing for Tracheostomy in Intubated Patients?

Bharat A. Panuganti, MD<sup>1</sup>; Philip Weissbrod, MD; Edward Damrose, MD<sup>2</sup>

Timing of Tracheostomy in Intensive Care Unit Patients

Terug/Vooruit

Tracheo-oesophageal fistula during  
tracheostomy and intermittent  
pressure ventilation<sup>1</sup>

S. HARLEY  
Hospital, Sully, Penarth, Glam.

Vorige/Volgende my  
the Critical

Charles C  
Jacqueli...

Tracheotomy: clinical review and guidelines<sup>\*</sup>

Paul De Leyn<sup>a,\*</sup>, Lieven Bedert<sup>b</sup>, Marion Delcroix<sup>c</sup>, Pieter Depuydt<sup>d</sup>,  
Geert Lauwers<sup>e</sup>, Youri Sokolov<sup>f</sup>, Alain Van Meerhaeghe<sup>g</sup>, Paul Van Schil<sup>h</sup>

<sup>a</sup> Department of Thoracic Surgery, University Hospital Leuven, Belgium  
<sup>b</sup> Department of Pneumology, Middelheim Hospital, Belgium  
<sup>c</sup> Department of Pneumology, University Hospital Leuven, Belgium  
<sup>d</sup> Department of Intensive Care, University Hospital Leuven, Belgium  
<sup>e</sup> Department of Thoracic and Vascular Surgery, ZOL Hospital Ghent, Belgium  
<sup>f</sup> Department of Thoracic Surgery, Erasme University Hospital Brussels, Belgium  
<sup>g</sup> Department of Pneumology, CHU A. Vesale, Montigny-le-Tilleul, Belgium  
<sup>h</sup> Department of Thoracic and Vascular Surgery, University Hospital Antwerp, Belgium

received in revised form 1 May 2007; accepted 24 May 2007; available online 27 June 2007

Tracheostomy Update  
When and How

Bradley D. Freeman, MD

To breathe or not to breathe: a review of artificial airway  
placement and related complications

John Donatelli · Ayushi Gupta · Ramya Santhosh · Todd R. Hazelton ·  
Leelakrishna Nallamshetty · Alvaro Macias · Carlos A. Rojas

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timing tracheostomy OR early AND late tracheostomy


Search

Advanced Create alert Create RSS

User Guide

RESULTS BY YEAR

TEXT AVAILABILITY

- Abstract
- Free full text
- Full text

ARTICLE ATTRIBUTE

- Associated data

ARTICLE TYPE

- Books and Documents
- Clinical Trial
- Meta-Analysis
- Randomized Controlled Trial
- 
- 

Filters applied: Meta-Analysis, Randomized Controlled Trial, Systematic Review, in the last 10 years, Adult: 19+ years. [Clear all](#)

Did you mean **timing of tracheostomy OR early AND late tracheostomy** (8 results)?

**Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review.**

Cite Adly A, Youssef TA, El-Beghermy MM, Younis HM.  
Eur Arch Otorhinolaryngol. 2018 Mar;275(3):679-690. doi:  
Share 10.1007/s00405-017-4838-7. Epub 2017 Dec 19.  
PMID: 29255970 Review.

**Incidence of sternal wound infection after tracheostomy in patients undergoing cardiac surgery: A systematic review and meta-analysis.**

Cite Toeg H, French D, Gilbert S, Rubens F.  
J Thorac Cardiovasc Surg. 2017 Jun;153(6):1394-7.  
Share 10.1016/j.jtcvs.2016.11.040. Epub 2016 Nov 17.  
PMID: 27964980 Free article

**Effect of early tracheostomy on mortality in patients with acute respiratory distress syndrome: a systematic review and meta-analysis.**

Cite ...  
Share ...  
PMID: ...

PUBLICATION DATE

- 1 year
- 5 years
- 10 years
- Custom Range

AGE

- Adult: 19+ years

Additional filters

**Early versus late percutaneous dilational tracheostomy in critically ill patients anticipated requiring prolonged mechanical ventilation.**

Cite Zheng Y, Sui F, Chen XK, Zhang GC, Wang XW, Zhao S, Song Y, Liu W, Xin X, Li WX.  
Share Chin Med J (Engl). 2012 Jun;125(11):1927-30.  
PMID: 22884055 Clinical Trial

**Early tracheostomy in patients with acute respiratory distress syndrome but has no impact on mortality: a retrospective study.**

Cite ...  
Share ...  
PMID: ...  
doi: ...  
Schmelzer T, Padberg ...

REVIEW ARTICLE

**Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review**  
Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Beghermy<sup>1</sup> · Hussein M. Younis<sup>1</sup>


**Early versus late tracheostomy in patients with acute respiratory distress syndrome: a randomized pilot trial SETPOINT.**

Cite ...  
Share ...  
PMID: 22264372 Clinical Trial.

**Early and late outcome after single step dilatational tracheostomy versus the guide wire dilating forceps technique: a prospective randomized clinical trial.**

Cite Fikkers BG, Staatsen M, van den Hoogen FJ, van der Hoeven JG.  
Share Intensive Care Med. 2011 Jul;37(7):1103-9. doi: 10.1007/s00134-011-2222-4. Epub 2011 Apr 12.  
PMID: 21484081 Free PMC article. Clinical Trial.

European Archives of Oto-Rhino-Laryngology (2018) 275:679–690  
<https://doi.org/10.1007/s00405-017-4838-7>

REVIEW ARTICLE

## Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Beghermy<sup>1</sup> · Hussein M. Younis<sup>1</sup>

### **Meta-analyse**

1987-2017

English

RCT

### **Groepen**

< 7 dagen

<14 dagen


<21 dagen

### **Medline search:**

timing of tracheostomy

early vs late tracheostomy

tracheostomy after intubation

European Archives of Oto-Rhino-Laryngology (2018) 275:679–690  
<https://doi.org/10.1007/s00405-017-4838-7>

REVIEW ARTICLE

Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

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**Table 1** Data collected from the articles involved in comparing early and late tracheostomy in adults

Author	Type of study	Patients with early tracheostomy					Patients with late tracheostomy						
		NO_	at (Days)	HAP	Mortality	MV Mean ± SD	ICU stay Mean ± SD	NO_	at (Days)	HAP	Mortality	MV Mean ± SD	ICU stay Mean ± SD
Ahmed and kuo [6]	Retrospective cohort	29	<7	11	4	15.7 ± 6	19 ± 7.7*	107	>7	14	1	20 ± 16	25.8 ± 11.8
Alali et al. [7]	Retrospective cohort	62	<8	–	48	10 ± 5.65*	13 ± 5.65*	59	>8	–	39	16 ± 6.36	19 ± 7
Alhajhusain et al. [8]	Retrospective cohort	39	<9	5*	5	15.1 ± 8.2*	16.6 ± 7.6*	63	>9	25	23	27.2 ± 10.9	27.2 ± 9.1
Arabi et al. [9]	Prospective cohort	29	<7	–	5	9.6 ± 1.2*	10.9 ± 1.2*	107	>7	–	15	18.7 ± 1.3	21 ± 1.3
Armstrong et al. [10]	Retrospective cohort	62	<7	–	7	–	15 ± 12*	95	>7	–	11	–	29 ± 26
Barquist et al. [11]	Prospective cohort	29	<8	28	2	8.57 ± 7.9	4.96 ± 6	31	>28	28	5	8.83 ± 9	5.26 ± 6.5
Ben-Avi et al. [12]	Retrospective cohort	90	<14	–	16*	16 ± 11	27 ± 21	109	>14	–	17	21 ± 14	31 ± 23
Bickenbach et al. [13]	Retrospective cohort	237	<10	–	53*	20.1 ± 11.4*	25.5 ± 15.7*	59	>10	–	24	29.5 ± 13.8	34.2 ± 24.1
Blot et al. [14]	Prospective cohort	61	<4	–	21	–	–	95	>14	–	20	–	–
Bosel et al. [15]	Prospective cohort	30	<3	–	4	15 ± 4.94	17 ± 6.36	30	7–14	–	2	12 ± 5.65	18 ± 7
Chen et al. [16]	Retrospective cohort	22	<21	–	–	44.4 ± 10.4*	–	38	>21	–	–	59.5 ± 9.6	–
Choi et al. [17]	Retrospective cohort	10	<10	4*	–	5.2 ± 6.5*	20.8 ± 6*	11	>10	9	–	29.2 ± 22.9	38 ± 18.5
Devarajan et al. [18]	Retrospective cohort	114	<10	63*	24*	–	526.5 ± 351.2*	114	14–28	71	46	–	698.7 ± 332.5
Flaatten et al. [19]	Retrospective cohort	230	<6	–	51*	–	–	231	>6	–	75	–	–
Ganuja et al. [20]	Retrospective cohort	101	<7	75	1	33.79 ± 29.62*	40.31 ± 26.49*	114	>7	83	4	40.73 ± 22.71	49.42 ± 19.24
Gatti et al. [21]	Prospective cohort	14	<7	–	–	12.6 ± 6.5*	15 ± 6.5*	10	>7	–	–	20.9 ± 10.7	27.1 ± 11.1
Gessler et al. [22]	Retrospective cohort	39	<7	19*	3*	17.4 ± 1.98*	–	109	>7	75	8	22.3 ± 1.77	–
Holloway et al. [23]	Retrospective cohort	24	<14	–	1	–	15 ± 7.5*	49	>14	–	1	–	19 ± 11.8
Hosseini et al. [24]	Retrospective cohort	13,386	<10	–	2837*	–	–	26,625	>10	–	7147	–	–
Huang et al. [25]	Retrospective cohort	11	<10	–	1*	–	16 ± 8*	27	>10	–	4	–	29 ± 13
Hyde et al. [26]	Retrospective cohort	53	<5	–	–	16.7 ± 11.4*	21.4 ± 11*	53	>5	–	–	21.9 ± 12.9	28.6 ± 16.3
Jeon et al. [27]	Retrospective cohort	39	<10	2*	2	11.4 ± 5.6*	19.9 ± 10.6*	86	>10	16	6	21.5 ± 15.5	31.1 ± 18.2
Keenan et al. [28]	Retrospective cohort	5402	<10	2479	568	–	–	4260	>10	2094	372*	–	–
Koch et al. [29]	Retrospective cohort	50	<4	19*	10	15.3 ± 5.4*	21.5 ± 7.5*	50	>4	32	11	21.1 ± 7.24	30.6 ± 7.5
Mahafza et al. [30]	Retrospective cohort	70	<21	20*	13*	–	–	26	>21	15	13	–	–
Mehta et al. [31]	Retrospective cohort	7595	<7	–	831	–	–	12,230	>7	–	1850	–	–
Moller et al. [32]	Retrospective cohort	81	<7	22*	–	12.2 ± 0.9*	16.7 ± 1*	104	>7	44	–	21.9 ± 1.3	26 ± 1.3
Pinheiro et al. [33]	Retrospective cohort	11	<8	6*	1*	–	–	17	>8	12	8	–	–
Puentes et al. [34]	Retrospective cohort	32	<7	–	8*	–	–	115	>7	–	34	–	–
Rizk et al. [35]	Prospective cohort	1577	<7	861*	–	–	–	1527	>7	1118	–	–	–
Rodrigues et al. [36]	Prospective cohort	51	<7	33*	–	12 ± 1*	16 ± 1*	55	>7	53	–	32 ± 3	37 ± 4
Rumbach et al. [37]	RCT	60	<2	5*	19*	7.6 ± 4*	4.8 ± 1.4*	60	>14	15	37	17.4 ± 5.3	16.2 ± 3.8
Scales et al. [38]	Retrospective cohort	1081	<9	–	426	–	–	4773	>13	–	1768	–	–
Tong et al. [40]	Case series	128	<7	3*	8*	21.47 ± 1.86*	17.52 ± 1.38*	464	>7	15	71	39.33 ± 1.33	26.27 ± 0.73
Villwock et al. [41]	Retrospective cohort	5591	<10	347*	699	–	29.1 ± 24.4*	7574	>10	644	894	–	36.8 ± 26
Villwock and Jones [42]	Retrospective cohort	53,749	<10	2429*	8947*	–	24 ± 9.6*	71,244	>10	3576	15,466	–	33 ± 11.6
Wang et al. [43]	Retrospective cohort	16	<10	7*	2*	14.9 ± 8.9*	–	50	>10	38	4	22.1 ± 7.6	–
Young et al. [44]	RCT	315	<4	–	139	–	13 ± 5.46	312	>10	–	141	–	13.1 ± 8.2
Yue et al. [45]	RCT	58	3	17*	–	–	–	61	15	30	–	–	–


# resultaten


**Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review**

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

# resultaten

## Incidence of HAP


Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

# resultaten

## Mortality



Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

# resultaten

## Duration of Mechanical Ventilaton



Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

# resultaten

## ICU Length of Stay

European Archives of Oto-Rhino-Laryngology (2018) 275:679–690  
<https://doi.org/10.1007/s00405-017-4838-7>

REVIEW ARTICLE

Timing of tracheostomy in patients with prolonged endotracheal intubation: a systematic review

Ahmed Adly<sup>1</sup> · Tamer Ali Youssef<sup>1</sup> · Marwa M. El-Begery<sup>1</sup> · Hussein M. Younis<sup>1</sup>

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







# Timing of Tracheostomy in ICU Patients: A Systematic Review and Meta-Analysis of Randomized Controlled Trials

Raffaele Merola <sup>1,\*</sup>, Carmine Iacovazzo <sup>1</sup>, Stefania Troise <sup>2</sup>, Annachiara Marra <sup>1</sup>, Antonella Formichella <sup>1</sup>, Giuseppe Servillo <sup>1</sup> and Maria Vargas <sup>1</sup>


## Bayesian analysis of a systematic review of early versus late tracheostomy in ICU patients

Laura Quinn<sup>1,2,\*</sup>, Tonny Veenith<sup>3</sup>, Julian Bion<sup>4</sup>, Karla Hemming<sup>1</sup>, Tony Whitehouse<sup>3</sup> and Richard Lilford<sup>1</sup>

<sup>1</sup>Institute of Applied Health Research, University of Birmingham, Birmingham, UK, <sup>2</sup>NIHR Birmingham Biomedical Research Centre, University Hospitals Birmingham NHS Foundation Trust and University of Birmingham, Birmingham, UK, <sup>3</sup>Department of Critical Care and Anaesthesia, Queen Elizabeth Hospital Birmingham, Birmingham, UK and

<sup>4</sup>Intensive Care Medicine, University of Birmingham, Birmingham, UK


## Tracheostomy: Epidemiology, Indications, Timing, Technique, and Outcomes

Nora H Cheung MD and Lena M Napolitano MD


In conclusion, there is no benefit to early tracheostomy in most ICU patients with acute respiratory failure, and waiting until  $\geq 10$  d of intubation and mechanical ventilation is recommended to determine whether ongoing respiratory support is required. Special patient populations

considered. Based on the available information, one could consider a tracheostomy as soon as it is apparent that weaning from artificial ventilation is unlikely to happen within two weeks after endotracheal intubation, in particular in neurological patients (Level C). Since this prediction is very difficult, tracheostomy generally should be delayed until at least 10 days after initiation<sup>52</sup>.

VS



## Percutaneous and surgical tracheostomy in critically ill adult patients: a meta-analysis


Christian Putensen<sup>1\*</sup>, Nils Theuerkauf<sup>1</sup>, Ulf Guenther<sup>1</sup>, Maria Vargas<sup>2</sup> and Paolo Pelosi<sup>2</sup>

### Key messages

- PT can be performed faster and reduces the risk for stoma inflammation and infection when compared with ST.
- PT is associated with increased technical difficulties when compared with ST.
- Among PT techniques, MDT and SSDT were associated with the lowest risk of intraprocedural technical difficulties and major intraprocedural bleeding.
- Risk for tracheal stenosis and odds for hospital survival were not different between PT techniques and ST.

# Chirurgische Tracheostomie



# Percutane tracheostomy



# Contra-Indicatie

absoluut

# Contra-Indicatie

relatief

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
- CUFF – SUBGLOTTISCH / SUBGLOTTIC

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# Toys for the boys



# Work-up



Contra-indicaties



Inspectie H/H



Echo hals



# Anatomie






Distribution of thyroid veins.

Distribution of thyroid arteries with associated laryngeal nerve, anterior view

F. Praxel 50

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
beademing

Dr. luchtweg

Dr. trach


glidescope

bronchoscoop

## Percutaneous dilatational tracheostomy versus fibre optic bronchoscopy-guided percutaneous dilatational tracheostomy in critically ill patients: a randomised controlled trial


Guanggui Shen<sup>1</sup> · Hongzhen Yin<sup>1</sup> · Yingya Cao<sup>1</sup> · Meijun Zhang<sup>1</sup> · Jingyi Wu<sup>1</sup> · Xiaogan Jiang<sup>1</sup> · Tao Yu<sup>1</sup> · Weihua Lu<sup>1</sup>


Received: 22 May 2018 / Accepted: 1 August 2018 / Published online: 11 August 2018  
© Royal Academy of Medicine in Ireland 2018

**Table 2** Incidence of complications


Characteristic	All ( <i>n</i> = 90)	PDT ( <i>n</i> = 45)	PDT-FOB ( <i>n</i> = 45)	<i>t</i> / $\chi^2$	<i>P</i>
Complication rate (%) ( <i>n</i> )	30% (27/90)	40% (18/45)	20% (9/45)	4.286	0.038
Minor haemorrhage controlled by compression (%) ( <i>n</i> )	26.7% (24/90)	33.3% (15/45)	20% (9/45)	2.045	0.153
Major haemorrhage requiring intervention (%) ( <i>n</i> )	1.1% (1/90)	4.4% (2/45)	0 (0/45)	–	0.494
Pneumothorax (%) ( <i>n</i> )	1.1% (1/90)	2.2% (1/45)	0 (0/45)	–	1.000
Total rate of first-time success (%) ( <i>n</i> )	78.9% (71/90)	64.4% (29/45)	93.3% (42/45)	11.27	0.001
Puncture with the puncture needle (%) ( <i>n</i> )	84.4% (76/90)	75.6% (34/45)	93.3% (42/45)	5.414	0.039
Insertion of the tracheostomy catheter (%) ( <i>n</i> )	94.4% (85/90)	88.9% (40/45)	100% (45/45)	5.294	0.056
Time of procedure (min)	11.3 ± 1.3	12.9 ± 1.1	9.8 ± 1.2	12.812	<0.001

*opereren*  
=  
*positioneren*

**TRACOE**<sup>®</sup>  
*twist*


**H+**  
**MC** Intensive Care





**Atos**  
Breathing-Speaking-Living



**TRACOE**<sup>®</sup>  
*twist*






**H+**  
MC Intensive Care

**I hate  
Complications**

**ORIGINAL ARTICLE**

# Safety and Feasibility of Very Early Bronchoscopy-assisted Percutaneous Dilatational Tracheostomy in Anterior Cervical Spine Fixation Patients


Amrutha Liz Paul<sup>1</sup>, Ram Varaham<sup>2</sup>, Kannan Balaraman<sup>3</sup>, S Rajasekaran<sup>4</sup>, Balasubramani VM<sup>5</sup>

Received on: 23 August 2022; Accepted on: 25 August 2022; Published on: 30 September 2022


**Table 4: Complications**



Bleeding (not amounting to transfusion)	4 (0.04%)
Bradycardia	7 (0.08%)
Desaturation	8 (0.09%)
Surgical-site infection	1 (0.01%)

Values presented as numbers (percentage)

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



## A

- Assess, Prevent, and Manage Pain

**EZ PICTURE BOARD BY VIDATAK**  
AN INNOVATION IN PATIENT COMMUNICATION

### ● I AM

short of breath	in pain	choking	feeling sick
hungry/thirsty	cold/hot	tired	dizzy
angry	afraid	frustrated	sad

### ● I WANT TO SEE

doctor	nurse	family	chaplain
--------	-------	--------	----------

### ● I WANT

to be suctioned	lip moistened	water	to be comforted	to sleep
tv/video/dvd	call light /remote	it quiet	lights off/on	to go home
to sit up	to lie down	to turn left/right	head of bed up/down	get out of bed

A	B	C	D	E	F	G	H	I	1	2	3	Thank You		
J	K	L	M	N	O	P	Q	R	4	5	6		I Love You	
S	T	U	V	W	X	Y	Z	.	7	8	9			❤️
'	,	?	!	SPACE							+	0	-	

SINGLE PATIENT USE. Please do not re-use between patients.

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
## Framework - English

<https://ilccare.org/the-fundamentals-of-care-framework/>

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## Patients' and nurses' experiences of fundamental nursing care: A systematic review and qualitative synthesis

Claire Pentecost BSc (Hons), PhD ✉, Julia Frost BA (Hons), MSc, PhD,  
Holly V. R. Sugg LLB Hons, MSc, PGDip, PhD, Angelique Hilli PhD ... See all authors ▾

Patient- and family-centered care interventions for improving the quality of health care: A review of systematic reviews


Myonghwa Park, Thi-Thanh-Tinh Giap ✉, Mihyun Lee, Hyun Jeong, Miri Jeong, Younghye Go

## Intensive care unit length of stay is reduced by protocolized family support intervention: a systematic review and meta-analysis

Systematic  
Volume 45

Contents lists available at ScienceDirect  
Intensive & Critical Care Nursing  
journal homepage: [www.sciencedirect.com/journal/intensive-and-critical-care-nursing](http://www.sciencedirect.com/journal/intensive-and-critical-care-nursing)

Review Article

### The effectiveness of family participation interventions for the prevention of delirium in intensive care units: A systematic review

Marli Lopo Vitorino <sup>a,b,\*</sup> ✉, Adriana Henriques <sup>a,b</sup>, Graça Melo <sup>a,b</sup>, Helga Rafael Henriques <sup>b</sup>

<sup>a</sup> University of Lisbon/Nursing School of Lisbon, Lisbon, Portugal

<sup>b</sup> Nursing Research, Innovation and Development Centre of Lisbon (CIDNUR), Nursing School of Lisbon, Lisbon, Portugal

Issues in Family

## An evaluation of family-centered care services a

Patient-centered care, nurse work environment and implicit rationing of nursing care in Swiss acute care hospitals: a cross-sectional multi-center study ☆

Stephanie Bachnick <sup>a</sup> ✉, Dietmar Ausserhofer <sup>a,b</sup> ✉, Marianne Baernholdt <sup>c</sup> ✉, Michael von <sup>a,d</sup> ✉, On behalf of the Match RN study group <sup>1</sup>

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# 04

## Decannulation ahead

From evidence to everyday practice in neurological patients  
Prof. Dr. Rainier Dziewas | Neuroloog Hospital Osnabrück  
President ESSD

A special thanks to ESSD

## ESSD 2026 16<sup>th</sup> Annual Congress

Advances in Measurement in Dysphagia:  
Innovation and Best Practices

The ESSD is delighted to host  
its annual congress in:

**Leiden, Netherlands**  
**From 19 to 23 October 2026**

<https://essd2026.org>

# Decannulation ahead: from evidence to everyday practice in neurological patients

---

Rainer Dzierwas

Department of Neurology and Neurorehabilitation

Klinikum Osnabrück - Academic Teaching Hospital of the University of Münster

# Outline

- > **General treatment options**
  - > **Guideline suggestions & basic principles & education**
- > The A<sup>2</sup>BC approach (I)
  - > Targeted diagnostic procedures
- > The A<sup>2</sup>BC approach (II)
  - > Targeted therapeutic interventions
- > The comprehensive decannulation algorithm
  - > Combining fast-track and standard-track pathways
- > Looking into the New England Journal of Medicine
  - > Fast-track decannulation works

# What do guidelines suggest?

Dziewas et al. *Neurological Research and Practice* (2021) 3:23  
<https://doi.org/10.1186/s42466-021-00122-3>

Neurological Research  
and Practice

**GUIDELINES**

**Open Access**

Diagnosis and treatment of neurogenic dysphagia – S1 guideline of the German Society of Neurology


# Looking into the guidelines

- > **Recommendation 22:** Patients with a tracheal cannula should be managed by a **multi-professional team**.
- > **Recommendation 23:** In tracheotomized patients with the therapeutic aim of decannulation, **swallowing function, oropharyngeal secretion management, vigilance and ability to cooperate, respiratory function, and airway anatomy, the voluntary and reflexive cough, as well as the amount, nature and clearing of the bronchial secretion,** should be evaluated regularly.
- > **Recommendation 24:** In tracheotomized patients, swallowing function should be evaluated with **FEES** and in particular the parameters **“secretion management”, “spontaneous swallowing rate”,** and **“laryngeal sensitivity”** should be investigated.


# Looking into the guidelines

- > **Recommendation 25:** In tracheotomized patients, the **location, fit and patency of the cannula, presence of granulation tissue and the placement of any existing fenestration** should be checked regularly.
- > **Recommendation 26:** If patients are intended to be weaned gradually from the tracheal cannula, a **physiological air flow through the upper airway** should be established to improve pharyngo-laryngeal sensitivity. Therefore, **the tracheal cannula's cuff should be intermittently deflated**, and the cannula be capped or a oneway speaking valve used.
- > **Recommendation 27:** During gradual weaning of the tracheal cannula, if necessary, in the clinical context, **the diameter of the inner cannula should be downsized** to reduce the airway resistance.
- > **Recommendation 28:** A definitive decannulation is usually possible if the cannula's cuff can be **continuously deflated with the cannula simultaneously being capped for 24-48 h** without complications.

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
Deut  
Ges  
[www.dgd.de](http://www.dgd.de)

# Outline

- > General treatment options
  - > Guideline suggestions & basic principles
- > The A<sup>2</sup>BC approach (I)
  - > Targeted diagnostic procedures
- > The A<sup>2</sup>BC approach (II)
  - > Targeted therapeutic interventions
- > The comprehensive decannulation algorithm
  - > Combining fast-track and standard-track pathways
- > Looking into the New England Journal of Medicine
  - > Fast-track decannulation works

# All in one: a recent review

Dziewas et al.

*Neurological Research and Practice* (2025) 7:18

<https://doi.org/10.1186/s42466-025-00376-1>

Neurological Research  
and Practice

REVIEW

Open Access

## Decannulation ahead: a comprehensive diagnostic and therapeutic framework for tracheotomized neurological patients

Rainer Dziewas<sup>1,3\*</sup> , Tobias Warnecke<sup>1</sup>, Bendix Labeit<sup>2</sup>, Volker Schulte<sup>1</sup>, Inga Claus<sup>3</sup>, Paul Muhle<sup>3</sup>, Anna Brake<sup>1</sup>, Lena Hollah<sup>1</sup>, Anne Jung<sup>3</sup>, Jonas von Itter<sup>3</sup> and Sonja Suntrup-Krüger<sup>3</sup>

# Targeted diagnostics: A<sup>2</sup>BC-Criteria

Category	Methods	Criteria
Airway safety	<ul style="list-style-type: none"><li>• FEES</li><li>• Secretion management</li><li>• Sensation</li><li>• Swallowing function</li></ul>	SESETD = 3 points
Airway anatomy	<ul style="list-style-type: none"><li>• Measurement of P<sub>IT</sub></li><li>• Endoscopy of upper and lower airway</li></ul>	P <sub>IT</sub> < 5 cm H <sub>2</sub> O No sign of critical airway obstruction
Bronchial secretions	<ul style="list-style-type: none"><li>• Suctioning frequency</li><li>• semiquantitative scores</li></ul>	Suctioning frequency ≤ 2/8h m-sqAS < 4 points
Cough strength	<ul style="list-style-type: none"><li>• Peak-flow measurement</li><li>• Semiquantitative scores</li></ul>	PCF ≥ 160 l/min, MEP ≥ 40 cmH <sub>2</sub> O SCSS ≥ 3 points; positive WCT

# SESETD-Algorithm

Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation

## Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation in Critically Ill Neurologic Patients

Tobias Warnecke, MD, PhD; Sonja Suntrup, MD; Inga K. Teismann, MD, PhD; Christina Hamacher; Stephan Oelenberg; Rainer Dziewas, MD, PhD

Warnecke et al. *Neurological Research and Practice* (2020) 2:9  
<https://doi.org/10.1186/s42466-020-00055-3>

Neurological Research  
and Practice

RESEARCH ARTICLE

Open Access

Inter-rater and test-retest reliability of the “standardized endoscopic swallowing evaluation for tracheostomy decannulation in critically ill neurologic patients”

Tobias Warnecke<sup>1†</sup>, Paul Muhle<sup>1,2\*†</sup>, Inga Claus<sup>1</sup>, Jens B. Schröder<sup>1</sup>, Bendix Labeit<sup>1,2</sup>, Sriramya Lapa<sup>3</sup>, Sonja Suntrup-Krueger<sup>1,2</sup> and Rainer Dziewas<sup>1</sup>

Muhle et al. *Neurological Research and Practice* (2021) 3:26  
<https://doi.org/10.1186/s42466-021-00124-1>

Neurological Research  
and Practice

RESEARCH ARTICLE

Open Access

Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation in Critically Ill Neurologic Patients – a prospective evaluation

Paul Muhle<sup>1,2\*</sup> , Sonja Suntrup-Krueger<sup>1,2</sup>, Karoline Burkhardt<sup>3</sup>, Sriramya Lapa<sup>4</sup>, Mao Ogawa<sup>5</sup>, Inga Claus<sup>1</sup>, Bendix Labeit<sup>1,2</sup>, Sigrid Ahring<sup>1</sup>, Stephan Oelenberg<sup>1</sup>, Tobias Warnecke<sup>1</sup> and Rainer Dziewas<sup>6</sup>

# SESETD-Algorithm


Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation

## Procedure:

1. First view

2. Suctioning & Cuff deflation & Speaking valve

3. SESETD protocol

# SESETD-Algorithm

Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation


- > Prospective observational study
- > FEES 4,75 d after weaning
- > Clinical swallow evaluation blinded to FEES
- > Main findings in FEES:
  - > Saliva pooling with aspiration (42)
  - > No spontaneous swallow (29)
  - > Severe pharyngeal hypesthesia (20)
- > Decannulation after FEES:
  - > 54 patients
  - > Necessity of recannulation: 1 patient
- > Decannulation deemed possible by CSE: 29 patients

**TABLE 1. Description of the Study Population**

Patient Characteristics	
Total, <i>n</i>	100
Sex: male/female, <i>n</i>	55/45
Age (yr)	56.4 ± 14.9
Diagnoses	
Cerebral infarction, <i>n</i>	51
Intracerebral hemorrhage, <i>n</i>	14
Guillain-Barre syndrome, <i>n</i>	11
Meningoencephalitis, <i>n</i>	6
Subarachnoid hemorrhage, <i>n</i>	3
Other, <i>n</i>	15
Duration of mechanical ventilation (d)	20.2 ± 13.1

# SESETD-Algorithm

Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation


> Decannulation:  
>  $227/377 = 60.2\%$


> Decannulation Failure:  
> Immediate:  $3/59 = 5.1\%$   
> Delayed:  $5/168 = 3.0\%$   
> Total:  $8/227 = 3.6\%$

# SESETD-Algorithm



Standardized Endoscopic Swallowing Evaluation for Tracheostomy Decannulation



# Patient case (1)


## Patient case (2)

# Patient case (2)

# Patient case (2)


# Outline

- > General treatment options
  - > Guideline suggestions & basic principles
- > The A<sup>2</sup>BC approach (I)
  - > Targeted diagnostic procedures
- > The A<sup>2</sup>BC approach (II)
  - > Targeted therapeutic interventions
- > The comprehensive decannulation algorithm
  - > Combining fast-track and standard-track pathways
- > Looking into the New England Journal of Medicine
  - > Fast-track decannulation works

# Treatment to target: A<sup>2</sup>BC-Criteria



Category	Methods
<b>Airway safety</b> <ul style="list-style-type: none"><li>• Secretion management</li><li>• Sensitivity</li> <li>• Swallowing function</li></ul>	<ul style="list-style-type: none"><li>• Saliva reduction (e.g. anticholinergics, botulinum toxin injection)</li><li>• Tactil-thermal, gustatory and olfactory stimulation</li><li>• Ice-chip protocol</li><li>• Above-cuff vocalization</li><li>• Pharyngeal electrical stimulation</li><li>• Oropharyngeal strengthening and coordination exercises</li><li>• Pharyngeal electrical stimulation</li></ul>
<b>Airway anatomy</b>	<ul style="list-style-type: none"><li>• Reduction of the TC outer diameter</li><li>• Customization of TC configuration</li><li>• Systemic/inhaled corticosteroids</li><li>• Minimally invasive or surgical interventions</li></ul>
<b>Bronchial secretions</b>	<ul style="list-style-type: none"><li>• Secretolysis (e.g. inhalation, secretolytics, fluid administration)</li><li>• Secretion reduction (e.g. anticholinergics)</li><li>• Secretion clearance (Air-stacking, oscillation therapy, mechanical cough aid, bronchoscopy)</li></ul>
<b>Cough strength</b>	<ul style="list-style-type: none"><li>• Breathing exercises (e.g. incentive spirometry)</li><li>• Cough reflex training (e.g. aerosolized capsaicin)</li></ul>

# Pharyngeal electrical stimulation (PES)


- > 3 consecutive days
- > 10 minutes of treatment per day

# Where it all started - PHAST-TRAC


Trial	PES	Sham		P
Suntrup et al.	15 (75.0)	2 (20.0)		0.0085
PHAST-TRAC	17 (50.0)	3 (8.8)		<0.001
<b>All</b>	<b>32 (59.3)</b>	<b>5 (11.4)</b>		<0.001


0.50 1.0 5.0 80.0  
Favours Sham Favours PES

# Providing real world evidence - PHADER

- Group A: Stroke, non-ventilated
- Group B: Stroke, ventilated

- Group C: Non-stroke, ventilated
- Group D: Traumatic Brain Injury

# PHADER - tracheotomized patients

Neurotherapeutics 21 (2024) e00433

Contents lists available at ScienceDirect

Neurotherapeutics

journal homepage: [www.sciencedirect.com/journal/neurotherapeutics](http://www.sciencedirect.com/journal/neurotherapeutics)

Original Article

Predictors of pharyngeal electrical stimulation treatment success in tracheotomised stroke patients with dysphagia: Secondary analysis from PHADER cohort study

Ivy Cheng<sup>a,b,c,g</sup>, Philip M. Bath<sup>d,e</sup>, Shaheen Hamdy<sup>b,f</sup>, Paul Muhle<sup>c,g</sup>, Satish Mistry<sup>f</sup>, Rainer Dziewas<sup>g,h</sup>, Sonja Suntrup-Krueger<sup>c,g,\*</sup>



Cheng et al. *Neurological Research and Practice* (2025) 7:23  
<https://doi.org/10.1186/s42466-025-00380-5>

Neurological Research and Practice

RESEARCH ARTICLE Open Access

Clinical predictors of outcome after pharyngeal electrical stimulation (PES) in non-stroke related neurogenic dysphagia after mechanical ventilation and tracheotomy: results from subgroup analysis of PHADER study

Ivy Cheng<sup>1,2,3,7</sup>, Philip M. Bath<sup>4,5</sup>, Shaheen Hamdy<sup>2,6</sup>, Paul Muhle<sup>3,7</sup>, Satish Mistry<sup>6</sup>, Rainer Dziewas<sup>7,8</sup> and Sonja Suntrup-Krüger<sup>3,7\*</sup>


# PHADER - Secondary analysis (1)

Baseline participant characteristics and intervention characteristics by group. Data are presented in mean (standard deviation), number (%) or median [interquartile range].

	All ventilated stroke patients (n = 98)	Patients who received PES during tracheotomy (n = 60)
<b>Participant characteristics</b>		
Age	66.6 (13.0)	66.6 (11.2)
Sex (Male / Female)	72 (73.5) / 26 (26.5)	44 (73.3) / 16 (26.7)
Body Mass index (BMI; kg/m <sup>2</sup> )		
by age groups		
50 years old or younger	24.6 (4.8)	25.8 (4.9)
51-60 years old	26.8 (9.1)	24.2 (4.6)
61-70 years old	26.2 (5.4)	26.1 (5.4)
71-80 years old	28.0 (4.4)	28.8 (3.6)
80 years old or older	26.5 (6.0)	26.4 (5.6)
Feeding status at baseline		
Oral, normal	0 (0.0)	0 (0.0)
Oral, supervision	1 (1.0)	0 (0.0)
Oral, with support	0 (0.0)	0 (0.0)
NGT or NJT	71 (72.4)	46 (76.7)
PEG or RIG	24 (24.5)	13 (21.7)
Other routes	2 (2.0)	1 (0.0)
NIHSS (/42)	13.4 (5.8)	13.5 (5.0)
Stroke type		
Ischaemic	74 (75.5)	46 (76.6)
Haemorrhagic	24 (24.5)	14 (23.3)
Stroke lesion location		
Supratentorial	84 (85.7)	54 (90.0)
Infratentorial	14 (14.3)	6 (10.0)
<b>Intervention characteristics</b>		
PES perceptual threshold at baseline (mA)	15.9 (7.9)	14.5 (7.6)
PES intensity on session 1 (mA)	30.9 (10.2)	30.0 (10.6)
Time from stroke onset to treatment (days)	29.5 [34.8]	32.5 [25.5]
<b>DSRS</b>		
Baseline	11.7 (1.2)	11.7 (1.2)
Day 5	10.8 (2.4)	11.1 (1.9)
Day 9	8.9 (3.8)	9.3 (3.6)
Day 92	5.3 (5.0)	5.6 (5.2)

DSRS: dysphagia severity rating scale [24]; NGT: nasogastric tube; NIHSS: National Institute Health Stroke Scale; NJT: nasojejunal tube; PEG: percutaneous endoscopic gastrostomy tube; PES: pharyngeal electrical stimulation; RIG: radiographically inserted gastrostomy tube.

# PHADER - Secondary analysis (2)

Participant characteristics	All ventilated non-stroke (n = 57)
Age	63.6 (15.5)
Sex (Male / Female)	40 (70.2) / 17 (29.8)
<i>Neurological conditions causing dysphagia</i>	
Traumatic brain injury	22 (38.6)
Critical illness polyneuropathy	15 (26.3)
Hypoxia	3 (5.3)
Seizures	3 (5.3)
Encephalitis	2 (3.5)
Guillain-Barré	2 (3.5)
Meningitis	2 (3.5)
Tumour	2 (3.5)
Brain abscess	1 (1.8)
Cavernoma	1 (1.8)
Cerebral oedema	1 (1.8)
Encephalopathy	1 (1.8)
Multiple sclerosis	1 (1.8)
Neurosarcoidosis	1 (1.8)
<i>Feeding status at baseline*</i>	
Oral, normal	0 (0.0)
Oral, supervision	0 (0.0)
Oral, with support	0 (0.0)
NGT or NJT	29 (50.9)
PEG or RIG	26 (45.6)
Other routes	2 (3.5)

**Table 2** Multiple linear regression findings for participants with neurogenic (non-stroke) dysphagia who required mechanical ventilation and tracheotomy (n = 57)

Participant characteristics	Change In DSRS at Day 5			Change In DSRS at Day 9			Change In DSRS at 3 months		
	$\beta$ [95% CI]	SE	p	$\beta$ [95% CI]	SE	p	$\beta$ [95% CI]	SE	p
Age	0.008 [-0.028, 0.045]	0.018	0.463	0.008 [-0.058, 0.074]	0.033	0.803	0.004 [-0.097, 0.105]	0.050	0.938
Sex	0.473 [-0.758, 1.704]	0.614	0.771	1.587 [-0.618, 3.792]	1.099	0.155	0.348 [-2.884, 3.580]	1.606	0.829
Intervention characteristics	Change In DSRS at Day 5			Change In DSRS at Day 9			Change In DSRS at 3 months		
	$\beta$ [95% CI]	SE	p	$\beta$ [95% CI]	SE	p	$\beta$ [95% CI]	SE	p
PES perceptual threshold at the first session	0.091 [-0.050, 0.232]	0.070	0.202	0.281 [0.034, 0.528]	0.123	<b>*0.027</b>	0.199 [-0.127, 0.526]	0.162	0.225
PES stimulation intensity at the first session	-0.012 [0.106, 0.083]	0.047	0.801	-0.096 [-0.263, 0.071]	0.083	0.254	0.000 [-0.218, 0.217]	0.108	0.998
Time from diagnosis to first PES	0.001 [-0.001, 0.003]	0.001	0.343	0.003 [0.000, 0.007]	0.002	0.077	0.013 [0.004, 0.022]	0.004	<b>*0.004</b>



CI: confidence interval; DSRS: dysphagia severity rating scale; PES: pharyngeal electrical stimulation; SE: standard error

> ...our findings suggested that in patients with neurogenic dysphagia, regardless of the lesion location, nature or chronicity, **early PES treatment can facilitate re-establishment of neural networks necessary to drive functional recovery.**


# Outline

- > General treatment options
  - > Guideline suggestions & basic principles
- > The A<sup>2</sup>BC approach (I)
  - > Targeted diagnostic procedures
- > The A<sup>2</sup>BC approach (II)
  - > Targeted therapeutic interventions
- > **The comprehensive decannulation algorithm**
  - > **Combining fast-track and standard-track pathways**
- > Looking into the New England Journal of Medicine
  - > Fast-track decannulation works

# The final challenge: Decannulation Algorithm

# All in one

# Outline

- > General treatment options
  - > Guideline suggestions & basic principles
- > The A<sup>2</sup>BC approach (I)
  - > Targeted diagnostic procedures
- > The A<sup>2</sup>BC approach (II)
  - > Targeted therapeutic interventions
- > The comprehensive decannulation algorithm
  - > Combining fast-track and standard-track pathways
- > Looking into the New England Journal of Medicine
  - > Fast-track decannulation works


# Does rapid

ESTABLISHED IN 1812

## High-Flow for

Gonzalo Hernández Martínez  
Ramón Ortiz, M.D.,  
Laura Colinas, M.D., Ph.D.  
Marina Garcia-de

# REDECAP Trial (1)

## Reducing decannulation times limiting capping

### > Inclusion criteria:

- > Patients with a tracheal cannula, weaned from ventilation

### > Exclusion criteria:

- > Unconsciousness (GCS <6)
- > Severe dysphagia
  - > 50 ml Water-Swallow-Test
  - > if pathological FEES according to SESETD protocol
- > Impaired airway patency
  - > Failed 5 minutes tube occlusion
  - > concordant findings in endoscopy
- > Severe neuromuscular disease (other than ICUAW)

# REDECAP Trial (2)

## Reducing decannulation times limiting capping

### > Study intervention:

#### > Intervention group:

- > High-flow oxygen and unblocked fenestrated tracheal cannula
- > Decannulation criterion: <2 suctioning episodes/8 hours for 24 hours

#### > Control group:

- > High-flow oxygen and unblocked fenestrated tracheal cannula & intermittent capping trials
- > Decannulation criterion: Successful 24 h capping trial


### > Primary endpoint:

- > Time to decannulation

### > Secondary endpoints:

- > Decannulation failure
- > Respiratory infections
- > Length of stay on the ICU
- > ICU readmissions
- > In-hospital death

# Primary Outcome

## No. of Patients

Intervention group	169	20	7	3	1
Control group	161	46	17	7	1

**Table 2. Primary and Secondary Outcomes.\***

**Prim**

Outcome	Control Group (N=161)	Intervention Group (N=169)	Difference (95% CI)
Primary outcome: median time to decannulation (IQR) — days†	13 (11 to 14)	6 (5 to 7)	7 (5 to 9)
Secondary outcomes			
Decannulation failure — no. (%)	9 (5.6)	4 (2.4)	3.2 (-1.2 to 8.1)
Weaning failure — no. (%)‡	27 (16.8)	11 (6.5)	10.3 (3.4 to 17.4)
Pneumonia — no. (%)	16 (9.9)	7 (4.1)	5.8 (0.2 to 11.8)
Tracheobronchitis — no. (%)	47 (29.2)	32 (18.9)	10.3 (1.0 to 19.3)
Median duration of stay (IQR) — days			
In the ICU§	35 (27 to 51)	32 (25 to 43)	3 (-1 to 11)
In the hospital	62 (38 to 105)	48 (33 to 71)	14 (9 to 33)
Death — no. (%)			
In the ICU	0	0	0 (-2.2 to 2.3)
In the hospital	8 (5.0)	4 (2.4)	2.6 (-1.7 to 7.4)
Sepsis — no. (%)	12 (7.5)	12 (7.1)	0.3 (-5.5 to 6.3)
Multiorgan failure — no. (%)	6 (3.7)	2 (1.2)	2.5 (-1.1 to 6.8)
Exploratory outcomes			
Decannulation before ICU discharge — no. (%)	104 (64.6)	139 (82.2)	-17.7 (-26.8 to -8.1)
Capping-trial failure — no. (%)¶	118 (73.3)	NA	NA
Median duration of stay (IQR) — days			
In the hospital after randomization	37 (20 to 66)	23 (14 to 36)	14 (10 to 31)
In the hospital after ICU discharge	27 (11 to 53)	16 (7 to 27)	11 (4 to 20)
ICU readmission — no. (%)	17 (10.6)	10 (5.9)	4.6 (-1.4 to 10.9)
Swallowing dysfunction at decannulation — no. (%)	16 (9.9)	15 (8.9)	1.1 (-5.4 to 7.6)

# Summary (1)

- > Tracheostomy is one of the most frequent procedures on the ICU and particularly common in severely affected patients.
- > Dysphagia is often the key reasons why decannulation has to be delayed in this scenario.
- > According to recent guideline recommendations, FEES should be used for bedside dysphagia evaluation in tracheotomized patients.
- > Restoring physiological airflow through the upper airway constitutes an important therapeutic principle during the rehabilitation of tracheotomized patients.
- > Targeted diagnostic and therapeutic interventions need to focus on the A<sup>2</sup>BC criteria
  - > Airway safety
  - > Airway patency
  - > Bronchial secretions
  - > Cough strength




## Summary (2)

- > The SESETD algorithm assesses airway safety and airway patency and is an early prognosticator for „time-to-decannulation“
- > Pharyngeal electrical stimulation (PES) shortens time to decannulation in stroke patients and other patient groups.
- > The comprehensive decannulation algorithms distinguishes fast-track and standard-track pathways.
- > A recent RCT showed that the option of fast-track decannulation reduces
  - > time-to-decannulation
  - > Complications
  - > length of stay parameters
  - > the decannulation failure rate.

Thank you!

[rainer.dziewas@klinikum-os.de](mailto:rainer.dziewas@klinikum-os.de)  
[dziewas@uni-muenster.de](mailto:dziewas@uni-muenster.de)

A special thanks to ESSD

## ESSD 2026 16<sup>th</sup> Annual Congress

Advances in Measurement in Dysphagia:  
Innovation and Best Practices

The ESSD is delighted to host  
its annual congress in:

**Leiden, Netherlands**  
**From 19 to 23 October 2026**

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# 05

## Beyond the voice

A holistic approach to laryngeal rehabilitation

Prof. Dr. Sarah Wallace | Professor & Clinical Chair NHS

Manchester University

# Above Cuff Vocalisation: Beyond the voice, a holistic approach to laryngeal rehabilitation

Professor Sarah Wallace OBE FRCSLT

Consultant Speech and Language Therapist

Wythenshawe Hospital, Manchester University NHS Trust

Honorary Clinical Chair, The University of Manchester

National Tracheostomy Safety Project

RCSLT Specialist Advisor

Visiting Professor - Manchester Metropolitan University & University of Leuven, Belgium

# Disclosures

FUNDED BY

**NIHR** | National Institute for  
Health and Care Research

1. ATOS Tracheostomy Advisory Board
2. Safe and Effective Above Cuff Tracheostomy Vocalisation (SEACTV) study funded by the NIHR (i4i Product Development Award)

Clinicaltrials.gov (NCT04647786) on 20/11/2020

Favourable ethical opinion on 4/3/2022 (IRAS-ID 278006; REC-Ref 22/NW/0022)

# Overview and aims

1

Effects of intubation and tracheostomy on integrated laryngeal function

2

Negative impact on patient well-being and outcome

3

ACV and laryngeal rehab

4


Recent ACV research

# Why laryngeal function matters

Larynx is a  
multi/shared organ  
system

Voice, swallowing, airway protection,  
cough, breathing

Function requires  
normal anatomy and  
physiology

Upper airway airflow,  
subglottic pressure,  
sensation and timing  
are vital


# Recovery milestones

Eating a proper meal

Decannulation

Taking the first  
sip of water

Talking to  
loved ones

# Aetiologies of laryngeal dysfunction

Intubation

Mechanical ventilation  
Respiratory disease

Weakness

Critical illness myopathy =  
91% dysphagic\*

Comorbidities

Neurology, respiratory, cancer,  
surgery, frailty

Delirium, cognition

Laryngo pharyngeal  
Reflux – laryngeal injury

# Intubation and tracheostomy stats

13-20 million people are intubated every year worldwide, ICU survival and post-airway sequelae rates are increasing

Laryngeal injury is associated with worse breathing and vocal symptoms

15-20,000 new tracheostomies per year in UK. 250,000 total per year in resource rich countries


Mortality rates in trachy pts - 10-60% dependent on comorbidity



Tracheostomy does not eliminate laryngeal risk, it adds its own set of complications

# Effects of intubation on laryngeal function

1. Laryngeal injury (in up to 83%), clinically significant in 10–20%
2. Desensitisation (in up to 77%), impairs cough, LAR
3. Disuse atrophy of laryngopharyngeal musculature
4. Sedatives associated with pharyngeal weakness and aspiration
5. Dysphagia and aspiration (in up to 60%)
6. Dysphonia (in up to 76%), a marker for laryngeal dysfunction
7. Reduced UES sphincter compliance
8. ↑ injury with duration, reintubation, tube size, cuff pressure, proning, sedatives

# Intubation related laryngeal injury

Mucosal ischemia

Vocal fold immobility - cuff pressure on RLN

Persistent oedema

Ulceration, granulomas

Arytenoid cartilaginous trauma

Posterior glottic stenosis - airway obstruction

BJA Education, 21(7): 250–257 (2021)  
doi: 10.1016/j.bjoe.2021.02.005  
Advance Access Publication Date: 21 April 2021

**Laryngeal complications after tracheal intubation and tracheostomy**

S. Wallace<sup>1</sup> and B.A. McGrath<sup>1,2,\*</sup>

<sup>1</sup>Wythenshawe Hospital, Manchester University NHS Foundation Trust, Manchester, UK and <sup>2</sup>University of Manchester Academic Critical Care, Wythenshawe Hospital, Manchester, UK

\*Corresponding author: [brendan.mcgrath@manchester.ac.uk](mailto:brendan.mcgrath@manchester.ac.uk)

**Keywords:** complications; dysphagia; dysphonia; intra-tracheal; intubation; larynx; speech and language therapy; tracheostomy

European Journal of Medical Research

Home About Articles Submission Guidelines [Submit manuscript](#)

Research | [Open access](#) | Published: 31 August 2024

**Incidence of post-extubation dysphagia among critical care patients undergoing orotracheal intubation: a systematic review and meta-analysis**

[Weixia Yu](#), [Limi Dan](#), [Jianzheng Cai](#), [Yuyu Wang](#), [Qingling Wang](#), [Yingying Zhang](#) & [Xin Wang](#)

*European Journal of Medical Research* 29, Article number: 444 (2024) | [Cite this article](#)

# Impact of intubation on laryngeal sensation

► Curr Otorhinolaryngol Rep. Author manuscript; available in PMC: 2014 Sep 22.  
Published in final edited form as: Curr Otorhinolaryngol Rep. 2013 Sep;1(3):178-182. doi: [10.1007/s40136-013-0018-5](https://doi.org/10.1007/s40136-013-0018-5)

**Neurophysiology and Clinical Implications of the Laryngeal Adductor Reflex**

[Amanda S Domer](#)<sup>1</sup>, [Maggie A Kuhn](#)<sup>2</sup>, [Peter C Belafsky](#)<sup>3</sup>,

Multicenter Study > Dysphagia. 2019 Aug;34(4):521-528. doi: [10.1007/s00455-019-09980-1](https://doi.org/10.1007/s00455-019-09980-1).  
Epub 2019 Jan 29.

**Relationship Between Laryngeal Sensation, Length of Intubation, and Aspiration in Patients with Acute Respiratory Failure**

[James C Borders](#)<sup>1</sup>, [Daniel Fink](#)<sup>2</sup>, [Joseph E Levitt](#)<sup>3</sup>, [Jeffrey McKeehan](#)<sup>4</sup>, [Edel McNally](#)<sup>5</sup>, [Alix Rubio](#)<sup>6</sup>

**Laryngeal Adductor Reflex**  
Bilateral firing of thyroarytenoid muscles to close VC & protect airway

LAR is absent in up to 77% post-extubation, LAR reduces with age

Diminished LAR is associated with intubation duration, oedema, secretions, dysphagia, aspiration

Oedema raises LAR threshold further

Severity of sensory impairment and the impact of airflow is not known

Aviv. *Ann Otol Rhinol Laryngol*. 1993 Oct;102(10):777-80, Labeit et al. *Neurogastroenterol Motil*. 2019 Nov;31(11):e13690.

# Effect of tracheostomy on laryngeal function

Loss of upper airway / laryngeal airflow

Loss of voice



Desensitisation, no afferent swallow trigger

Laryngopharyngeal deconditioning & weakness

Breath / swallow pattern incoordination


# Effect of cuff inflation ‘use it or lose it’

- Loss of sensation - silent aspiration
- Loss of subglottic pressure → ineffective cough, pharyngeal residue
- Secretions pool in upper airway


# The swallow sensory-motor loop



# Breath– swallow coordination

“I can’t swallow right now, I have to keep breathing”

- Central Pattern Generators for respiration and swallowing are intertwined in the brainstem
- Sophisticated degree of coordination between the two systems
- Chewing disrupts respiratory cycles
- Larger lung volumes needed for liquids
- Hypercapnia reduces swallow frequency
- Ventilator alters synchrony with phonation, swallow apnoea

# Impact of tracheostomy on secretion management


**Increased secretion volume and viscosity**  
**Risks – prolonged weaning, ICU stay**

# Consequences of laryngeal dysfunction

- Prolongs wean & need for trache
- ↑ ICU/hospital LoS
- Malnutrition, NBM, thirst
- Aspiration pneumonia
- Chronic dysphagia, tube feeding
- Chronic dysphonia
- Worse psychological well-being
- ↑ mortality

# Psychological impact of voice loss

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ELSEVIER

What matters most to adults with a tracheostomy in ICU and the implications for clinical practice: a qualitative systematic review and metasynthesis

Helen Newman, MSc<sup>a,b,\*</sup>, Gemma Clunie, PhD<sup>c,d</sup>, Sarah Wallace, MSc<sup>e,f</sup>, Christina Smith, PhD<sup>g</sup>, Daniel Martin, PhD<sup>a,h</sup>, Natalie Pattison, PhD<sup>i,j</sup>

Check for updates

**Voice = IDENTITY**

Seen and heard as a unique person

A black and white close-up portrait of Morgan Freeman. He is looking slightly to the right of the camera with a subtle, thoughtful expression. The lighting is soft, highlighting the texture of his skin and the details of his facial features. The background is a plain, light gray.

**You are now reading this**

**in my voice**

# Voice loss impacts care


Difficulty optimising pain relief

3x more likely to have an adverse medical event

Difficulty managing delirium

Failed attempts to communicate induce staff stress


# Prevention and assessment of laryngeal dysfunction

## Prevention

- Smaller ETTs, careful cuff pressure monitoring
- Early tracheostomy
- Early laryngeal rehabilitation

## Assessment

- Clinical bedside exam
- FEES
- Detect injury, oedema, LPR, secretion issues
- Guide laryngeal rehabilitation

# Aims of laryngeal rehabilitation

1

Restore  
upper airway  
physiology

2

Improve  
laryngeal  
function and  
outcomes

3

Safer, faster  
weaning and  
decannulation

4

Individualised  
rehab and  
humanisation

5

Improve QoL

# Laryngeal rehab starts with restoring airflow

1. Early cuff deflation + One Way Valve (within 24-48hrs if haemodynamically stable)

2. ACV if cuff inflated

Check for updates

Journal of the Intensive Care Society

Journal of the Intensive Care Society  
2022, Vol. 0(0) 1-8  
© The Intensive Care Society 2022  
Article reuse guidelines:  
[sagepub.com/journals-permissions](https://sagepub.com/journals-permissions)  
DOI: 10.1177/17511437221113162  
[journals.sagepub.com/home/jics](https://journals.sagepub.com/home/jics)

SAGE

Special Article

**Benefits and options for voice restoration in mechanically ventilated intensive care unit patients with a tracheostomy**

Sarah Wallace<sup>1,2</sup>, Sue McGowan<sup>3</sup> and Anna-Liisa Sutt<sup>4,5,6</sup>


# What is Above Cuff Vocalisation (ACV)?

Retrograde airflow delivered above the cuff via subglottic port, when cuff is inflated, no ventilator adjustments needed

## **Rationale**

- Voice when cuff inflated
- Early laryngeal stimulation
- Restore physiology and sensation
- Precursor or adjunct to cuff deflation + OWV

# ACV Above Cuff Vocalisation

# From 'speaking' to laryngeal rehabilitation


ACV is about laryngeal rehabilitation, not just voice

Laryngeal dysfunction is a multidisciplinary concern

Successful weaning depends on the larynx

# ACV can support early laryngeal rehab

Provokes vocal fold closure, LAR

Stimulates swallow

Expels secretions

2-5l airflow/min

# Safety, voice and comfort

**Table 2.** Effectiveness and complications of ACV.

Original article

## Safety and feasibility of above cuff vocalisation for ventilator-dependant patients with tracheostomies

Brendan A McGrath<sup>1</sup> , Sarah Wallace<sup>2</sup>, Mark Wilson<sup>2</sup>, Leanne Nicholson<sup>2</sup>, Tim Felton<sup>1</sup>, Christine Bowyer<sup>1</sup> and Andrew M Bentley<sup>1</sup>

Nausea	1	1.1
Patient asked to remove	1	1.1
Total	91	100.0

ACV: above cuff vocalisation.

Fiberoptic endoscopic evaluation of swallowing (FEES)

27.5%

# Patient selection for ACV

1. Patent upper airway
2. Stable ventilation
3. Inflated cuff all or most of the time
4. Subglottic suction tracheostomy tube

## Contraindications

- Airway obstruction
- Surgical emphysema
- Intolerance

[Acta Anaesthesiol Scand](#). 2021 Jan; 65(1): 15–25.

PMCID: PMC7756796

Published online 2020 Nov 1. doi: [10.1111/aas.13706](https://doi.org/10.1111/aas.13706)

PMID: [32920849](https://pubmed.ncbi.nlm.nih.gov/32920849/)

## Above cuff vocalisation (ACV): A scoping review

[Antonija Petosic](#),<sup>1, 2</sup> [Marit F. Viravong](#),<sup>3</sup> [Anna M. Martin](#),<sup>4</sup> [Cecilie B. Nilsen](#),<sup>1</sup> [Kjell Olafsen](#),<sup>5</sup> and [Helene Berntzen](#)<sup>1</sup>

- 17 studies
- Audible voice or whisper in 88% across studies
- QoL scores improved (V-RQOL and QOL-MV)
- Minor complications in 27%
- Two serious adverse events: subcutaneous emphysema

Observational Study > Arch Phys Med Rehabil. 2022 Mar;103(3):394-401.

doi: 10.1016/j.apmr.2021.08.016. Epub 2021 Sep 22.

## Determining the Prevalence, Implementation Approaches, and Opinions of Above Cuff Vocalization: A Survey of Health Care Professionals

Claire S Mills<sup>1</sup>, Emilia Michou<sup>2</sup>, Mark C Bellamy<sup>3</sup>, Heidi J Siddle<sup>4</sup>, Cathy A Brennan<sup>5</sup>, Chris Bojke<sup>5</sup>

ACV practice varies widely  
Barriers - lack of knowledgeable staff,  
access to training

## Above Cuff Vocalization (ACV) Safety Guidelines

Above cuff vocalization (ACV) involves applying an external airflow through the subglottic port of a tracheostomy tube to restore laryngo-pharyngeal airflow. This can support vocalization and improve airway protection and swallowing.<sup>1</sup> However, clinical practice varies, safety concerns exist, and no clear safety guidelines are available.

A nominal group technique approach was used to achieve consensus on a range of safety statements with a group of eight international clinical and academic ACV experts. Participants included: one occupational therapist, one physician, one physiotherapist, and four speech and language therapists. Six participants were from the UK and one from Denmark. This study was funded and supported by Atos Medical, UK.

- 1 The professional(s) leading ACV assessment should be competent in complex dysphagia assessment, tracheostomy management, and assessment of voice and communication.
- 2 All staff assessing and delivering ACV should be trained and have completed ACV competencies.
- 3 Separate competencies are needed for assessment of ACV and delivery of ACV.
- 4 The contraindications for ACV are: complete lack of upper airway patency; bleeding from stoma; surgical emphysema.
- 5 The precautions for ACV are: pain on application; glottic/supraglottic bleeding; oedema in the pharynx/larynx; secretions emerging from the stoma site; partial lack of upper airway patency; tracheostomy tube not in optimal position; granulation tissue; previous history of surgical emphysema.
- 6 Assessment of airway patency and safety should involve: review of airway history; FEES/FNE if indicated; assessing patient comfort during ACV; monitoring for audible/tactile airflow on ACV application; gradual increase of flow rate starting at 1 L/min.
- 7 The upper limit for airflow delivery should be 5 L/min.
- 8 The maximum duration of airflow delivery should be 15 minutes per session.
- 9 Airflow delivery for ACV should only be delivered intermittently using a thumb port ensuring that airflow is stopped during
- 10 Patients should be continuously observed throughout ACV sessions.

ATOS ACV safety guidelines

# Modelling upper airway airflow on ACV


2.17 l/min matched natural exhalation velocity/pressures - enough for voice

Flow rates  $\geq 5$  l/min generate significant turbulence & pressure, no vocal benefit


# Latest ACV research

Manchester, device development

**Above Cuff  
Vocalisation**  
(ACV)

‘Safety and Effectiveness of  
Above Cuff Tracheostomy  
Vocalisation’

# SEACTV study aims

Test ACV safety parameters and measure effectiveness

Determine impact of airflow on voice, swallow, sensation, secretion management & aspiration severity

# Methods

Design and build prototype, bench test

Clinical trial on ICU

In-built data chip captures episodes, alerts

E-Stim catheter measures pharyngeal sensory thresholds  
Pre and post ACV


FEES  
Pre ACV, with ACV, and follow-up

Daily patient & staff ratings of voice, comfort & anxiety

E-Stim catheter

Visual Analogue Scale scores  
0 = no voice, uncomfortable  
100 = normal voice, comfortable

# ICU recruitment


n=59 pts

39 had SEACTV

20 had standard care  
(cuff up/ down, OWV, basic ACV)

Excluded cuff down, ACV unsafe

# Results- Safety and comfort

SEACTV device use - 5.3 days per patient, 37.5 minutes each episode

## Voice

35/39 (89.7%) achieved voice everyday, voice quality was rated 22-95/100

## Tolerance

Humidification increased the duration of ACV by 23 mins over 7 days



## Ease of use

72/100, highest rating 100

SEACTV delivered ACV for longer with no adverse effects  
SEACTV was effective, comfortable, and easy to use


# Results Patient anxiety

All patients had high anxiety levels, standard care pts stayed highly anxious  
SEACTV - anxiety levels reduced quickly, stayed low


# Results Laryngopharyngeal sensation

4 sensory thresholds recorded

All thresholds dropped = sensation improved  
Quantifies the impact of airflow on sensory impairment

# Results Aspiration severity

Immediate and significant improvement in aspiration severity with airflow  
Earlier oral intake

# FEES


Cuff inflated

Larynx full of secretions

ACV airflow starts, see  
bubbles  
Cough and swallow trigger  
spontaneously

Airway clears  
Vocal cords moving  
Laryngeal anatomy visible  
No secretions /aspiration  
Cough/swallow frequency  
improve

Switching the larynx back on with AC

# Take home messages

Laryngeal dysfunction is common and complex

ACV provides holistic laryngeal rehabilitation earlier

Supports humanisation and QoL

Research gaps – larger studies, protocol refinement

# Thank you, questions?

sarah.wallace@mft.nhs.uk

# Communication with ACV

Jeanette, 56yrs old Heart transplant

Before ACV and after ACV  
airflow is turned on

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