

Freevent Dressing Foam



Product description:

The Freevent Dressing Foam 65x70 is a single use foam dressing that provides a protection between the tracheostomy tube and the skin and absorb secretions. It is intended for patients with small secretion.

Document ID:	PF093-01-TechInfo	Edition:	3.0
Manufacturer:	Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden		
Classification: (EU) 2017/745	Class I, Rule 1		
Intended Use:	Freevent Dressing Foams are single use tracheostomy dressings that provide protection between the tracheal cannula and the skin and absorb secretions.		
Use specifications:	Intended medical indication Patients using a tracheostomy tube. Intended patient population Adult and pediatric patients using a tracheostomy tube. Intended to be used by medical personnel or patients with sufficient cognitive ability and manual dexterity who are judged as able to manage the device independently by a clinician. Intended usage Single use. Intended part of the body/type of tissue applied to or interacted with Skin around tracheostoma. Intended user profile Medical personnel or patients with sufficient cognitive ability and manual dexterity who are judged as able to manage the device independently by a clinician. Intended conditions of use FreeVent Dressings are for single use and for daily continuous use. Intended for use at home (indoor and outdoor), care facilities and hospitals, during rest, sleep and exercising under normal daily environment without any hygienic or environmental restrictions regarding temperature, moisture etc. Freevent Dressing Foam should be exchanged when needed, or at least every 24 hours.		
Operating principles	The FreeVent Dressing Foam is to be placed between the tracheal cannula and the skin. For hygienic reasons, the FreeVent Dressing Foam should be exchanged when needed/saturated (indicated by increase in size of the dressing foam), at least with every cannula change, or every 24 hours. The FreeVent Dressing Foam can either be used alone or in combination with other tracheostomy dressings.		
Contraindications:	None		
CE Mark:	Yes. Devices are CE-marked.		
GMDN code:	15624 (Tracheostomy tube dressing)		
Sterilization:	Non-sterile		

Product Information

Raw material:	Polyurethane foam
Latex information:	Not manufactured with natural rubber latex
Biological origin:	The device is not manufactured with materials derived from human or animal source.
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2 °C – 42 °C.
Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.
Hazardous components:	None
Expiration date:	5 years after manufacturing.
Packaging:	10 products are packed in a Mini-grip plastic bag made of polyethylene together with instructions for use.

Devices under Basic UDI-DI: 7331791-COM-0-000-0004-5B

REF	Name	UDI-DI
1431	Freevent Dressing Foam 65x70	07331791006784

Atos Medical AB compatible products:

None

Document Approvals
Approved Date: 2024-11-04

Task: Approval Task Verdict: Approve	HUABOD Adrienn Bodonyi, Product Support Specialist (huabod@coloplast.com) Issuer 24-Sep-2024 09:18:42 GMT+0000
Task: Approval Task Verdict: Approve	ADEL.KHWATMI Adel Khatmi, Sustaining Engineer (adel.khatmi-atosmedical@coloplast.com) Technical / Specialist 03-Nov-2024 16:19:19 GMT+0000
Task: Final Approval Verdict: Approve	KARNIL Karolina Nilsson, Head of Regulatory Affairs (karolina.nilsson-atosmedical@coloplast.com) Regulatory 04-Nov-2024 20:52:01 GMT+0000