# FREEVENT<sup>®</sup>

## Dressing Softfoam

**MD**



**CE** 0413

**STERILE R**

## Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Patents and trademarks

Freevent® is a registered trademark owned by Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

	Length (mm)	Width (mm)	Thickness (mm)
8034 Freevent Dressing Softfoam L	80	80	5
8035 Freevent Dressing Softfoam S	50	80	5
8036 FreeventDressing Softfoam Slim L	80	80	2,5

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Producent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Производитель; Produzen; Pengilang; 製造業者; 제조사; 製造商; 制造商; תוצר; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Üretim tarihi; წარმოების თარიღი; Ἀρτοποιημάτων ἀνάσπασις; İstehsal tarixi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조일자; 製造日期; 制造日期; תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Anvendes inden; Brukes innen; Viimeinen käyttöpäivä; Síðasti notkunardagur; Kölbl-

lik kuni; Derīguma termiņš; Naudoti iki; Datum expirace; Felhasználhatóság dátuma; Dátum najneskoršej spotreby; Rok uporabnosti; Data ważności; Data expirării; Rok uporabe; Rok upotrebe; Ημερομηνία λήξης; Срок на годност; Son kullanım tarihi; Յարցօրօս; Ἡμερομηνία λήξης; Son istifadə tarixi; Срок годности; Gunakan sebelum tanggal; Tarikh guna sebelum; 使用期限; 사용 기한; 使用截止日期; 使用截止日期;

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## Intended use

Freevent Dressings are sterile single use dressings that absorb secretion and provide protection for the skin around the tracheostomy tube for patients with a tracheostoma.

## Warnings

- When changing the dressing be cautious in order to prevent accidental de-cannulation.

## Precautions

- Freevent Dressings may be used on irritated skin. However, if skin irritation persists for more than 14 days, stop using the dressing and consult your clinician.
- The dressing should be changed at least every 24 hours or more often if needed to prevent skin irritation.
- Do not use the dressing if the pouch is damaged.
- Do not modify the dressing by eg cutting.
- Applying the wrong side of the dressing towards the skin may cause skin irritation.
- The dressings may irritate the skin. Stop using the dressing if skin irritation develops and consult your clinician.

## Description of the device

Freevent Dressing Softfoam is made of soft polyurethane foam that provides medium to high absorption of secretion.

## Instructions for use

### Operating instructions

Freevent Dressing Softfoam has a pre-perforated slit that can be torn open, as showed in Figure 1. Apply the dressing with the foam side against the skin and shiny side up. The Dressing should be placed between the tracheostomy tube and the skin with the slit upwards, see Figure 2.

When using a dressing without opening the slit, the tracheostomy tubes are inserted through the opening of the dressing. The tracheostomy tube is then inserted in the tracheostoma and fixed with a tube holder.

### Disposal

Freevent Dressings shall be changed at least every 24 h. Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

### Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Atos Medical AB, Sweden.