
PROVOX[®] 2

Instructions for Use
Clinician's

Rx
ONLY

CE 0413

STERILE EO

Atos

Atos Medical Your voice

Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7




Figure 8


Figure 9


Figure 10


Figure 11


Figure 12


Figure 13



Figure 14

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox is a registered trademark of Atos Medical AB. Provox2 is protected by US patent no 5.976.151, JP 4.122.058 and other patents.

Unregistered copy

Contents

ENGLISH	8
DEUTSCH	15
NEDERLANDS.....	23
FRANÇAIS	31
ITALIANO.....	39
ESPAÑOL	47
PORTUGUÊS.....	55
SVENSKA	63
DANSK.....	70
NORSK.....	77
SUOMI	84
ČESKY	91
POLSKI.....	98
HRVATSKI	106
ΕΛΛΗΝΙΚΑ	113
TÜRKÇE	122
РУССКИЙ	129
עברית.....	143
العربية	150
ORDERING INFORMATION...	151

Contents

1. Descriptive information.....	9
1.1 Intended use	9
1.2 CONTRAINDICATIONS.....	9
1.3 Description of the device	9
1.4 WARNINGS	9
1.5 PRECAUTIONS.....	10
2. Instructions for use.....	10
2.1 Voice prosthesis replacement	10
Preparation	10
Loading the voice prosthesis into the insertion tool.....	10
Inserting the new prosthesis into the puncture (Anterograde replacement)	10
Replacement by Overshooting	11
Retrograde replacement.....	11
2.2 Device lifetime.....	11
2.3 Prosthesis maintenance	12
2.4 Accessories	12
3. Adverse Events/Trouble Shooting Information.....	12
4. Date of Printing and User Assistance Information.....	14

1. Descriptive Information

1.1 Intended use

The Provox2 Voice Rehabilitation System is intended for use in surgical, prosthetic voice restoration after total laryngectomy. The prosthesis may be inserted by the physician at the time of the total laryngectomy (primary puncture), or at a later date (secondary puncture), or may be used to replace the present prosthesis.

1.2 CONTRAINDICATIONS

There are no known contraindications for use or replacement of the Provox2 prosthesis among patients already using prosthetic voice rehabilitation.

1.3 Description of the device

The Provox2 Voice Prosthesis is a sterile single use one-way valve that keeps a TE-puncture open for speech, while reducing the risk of fluids and food entering the trachea. It is made of medical grade silicone rubber and radio-opaque fluoroplastic. The prosthesis is not a permanent implant and needs periodic replacement. The outer diameter is 7.5 mm (22.5Fr). The Provox2 is available in lengths 4.5, 6, 8, 10, 12.5 and 15 mm.

The Provox2 Voice Prosthesis package contains:

- 1 Provox2 Voice Prosthesis, an Inserter and pre-lubricated Loading Tube (Fig. 1).
- 1 Provox Brush of size corresponding to the voice prosthesis, non-sterile
- 2 Provox2 Instructions for Use (Clinician's and Patient's)
- 1 Provox Brush Instructions for Use

Compatibility with MRI, X-ray and radiation therapy

All Provox voice prostheses except Provox ActiValve are MRI-safe (tested up to 3 Tesla), and safe to use during X-ray and radiation therapy (tested up to 70 Gy).

1.4 WARNINGS

Re-use and re-processing may cause cross-contamination and damage to the device, which could cause patient harm.

Dislodgement or extrusion of the Provox2 voice prosthesis from the TE puncture and subsequent ingestion, aspiration or tissue damage may occur. (See also Section 3).

To reduce the risk of dislodgment/extrusion:

- Select proper prosthesis size.
- If laryngectomy tubes or stoma buttons are used, choose devices that do not exert pressure on the prosthesis or catch onto the tracheal flange during insertion and removal.
- Instruct the patient to consult a physician immediately if there are any signs of tissue edema and /or inflammation /infection.
- Inform the patient to use only original Provox accessories of corresponding size for maintenance and to avoid all other kinds of manipulation of the voice prosthesis.

1.5 PRECAUTIONS

- During loading of the prosthesis make sure that the esophageal flange is completely folded forward in the loading tube (Fig. 7). Incorrect loading can cause device damage with immediate valve failure and fluid aspiration.
- Carefully assess any patient on anticoagulant treatment or with bleeding disorders for the risk of bleeding prior to prosthesis placement or replacement.

2. Instructions for use

2.1 Voice prosthesis replacement

Preparation

Determine the required length of the voice prosthesis by using the current (old) prosthesis as measuring device. If there is more than 3 mm / ~0.12 inches space between the tracheal flange and the mucosal wall, a shorter prosthesis should be used. If the prosthesis fit is too tight, a longer prosthesis should be used. Alternatively, the Provox Measure (see Accessories) can be used to determine the length.

Loading the voice prosthesis into the insertion tool:

1. Secure the Safety Strap of the prosthesis by threading it into the slit (Fig. 2) of the Inserter Pin. Fixate the safety strap by pulling it into the narrower part of the slit. (Fig. 3). **CAUTION:** Always attach the safety strap of the prosthesis to the Inserter Pin. It could otherwise be lost during insertion causing aspiration or ingestion.
2. Place the prosthesis on top of the Inserter with the tracheal side down (Fig.4).
3. Fold the esophageal flange forward by pinching it together between thumb and index finger. Hold the loading tube with your other hand. Insert the folded esophageal flange of the prosthesis into the slit of the loading tube (Fig. 5).
4. Push the Inserter Pin forward about 1 cm (~0.4 inches) (Fig.6).
5. Press the protruding part of the flange forward with your thumb into the loading tube. The esophageal flange should now be completely folded forward (Figure 7).
6. Push the inserter into the Loading Tube until mark line 1 is level with the edge of the loading tube (Fig.8). The Provox2 voice prosthesis is now prepared for placement.

CAUTION: If the esophageal flange is not completely folded forward, too much force will be needed to press the prosthesis forward which could cause device damage.

Inserting the new prosthesis into the puncture (Anterograde replacement).

Caution: Do not “inject” the prosthesis using the inserter like a syringe.

1. **Pull out the current prosthesis** with a non-toothed hemostat.
2. **Gently insert the Loading Tube** in the TE-puncture until the tip is completely inside the puncture.
3. **Insert the prosthesis** by holding the loading tube stable with one hand and push the inserter forward with your other hand until mark line 2 is level with the rim of the loading tube (Fig.9).

4. **Check the prosthesis placement** by holding the loading tube and inserter together and pulling them back through the puncture until the tracheal flange reaches the posterior esophageal wall (Fig. 10).
5. **Release the prosthesis** by removing the loading tube (the inserter remains in position) (Fig. 11). If the tracheal flange does not unfold automatically, rotate and pull the inserter pin, or use a non-toothed hemostat to pull the tracheal flange and rotate this in position.
6. **Check proper position** by rotating and lightly pulling the tracheal flange. Thereafter place the voice prosthesis so that the safety strap is pointing downwards.
7. **Check valve function** by asking the patient to speak (Fig. 12) and by observing the prosthesis while the patient drinks water (Fig. 13).
8. **Cut off the safety strap** using scissors on the outer edge of the tracheal flange (Fig. 14). The Provox2 voice prosthesis is now ready for use.

Note: If the insertion is not successful, it can be repeated with the same tools and devices as long as the safety strap has not been cut off. **If the procedure has to be repeated more than twice**, add additional water-soluble lubrication inside the loading tube.

Replacement by Overshooting

Overshooting means that during the insertion of the voice prosthesis it is initially pushed entirely into the esophagus. This approach can be useful when using very short prostheses (4.5 mm), if there is a very thin puncture wall or a mucosal fold in the esophageal front wall entraps the esophageal flange.

Procedure: Prepare and insert the valve as described in steps 1-3 above but insert the complete prosthesis into the esophagus, by pushing the Inserter Pin past mark line 2. Release the prosthesis as described in step 5. Leave the safety strap attached to the inserter.

Use a non-toothed hemostat to pull the prosthesis into position. **CAUTION:** Do not pull by the safety strap. Slight rotation might be helpful to fully position the prosthesis. The end of the tracheal flange should point downwards into the trachea and the identification knobs must always face the stoma. Cut the safety strap only after the correct fit is obtained.

Retrograde replacement

Provox GuideWire (REF 7215) may be used for retrograde insertion of the voice prosthesis in cases where the Provox2 cannot be inserted in an anterograde manner. See Instructions for Use accompanying the Provox GuideWire. Retrograde insertion must be performed by a clinician.

2.2 Device lifetime

Biofilm (bacteria and candida) overgrowth will occur in almost all patients. This will eventually cause material degradation, leakage through the prosthesis or other incompetence of the valve. Radiotherapy, salivary content, gastroesophageal reflux, and dietary habits may influence how rapidly biofilm forms. It is therefore not possible to predict the device lifetime.

Laboratory testing shows that in the absence of bacteria and yeasts the device will maintain its structural integrity for 12 months of simulated use. Usage beyond this limit is under the sole discretion of the prescriber.

Indications for replacement include leakage through or (partial) blockage of the valve, degradation of the materials, or medical indications such as problems with the puncture tract (See also Section 3).

2.3 Prosthesis maintenance

Instruct the patient to clean the prosthesis at least twice a day and after each meal.

Use the **Provox Brush** and if possible also the **Provox Flush** to clean the prosthesis from food remnants and mucus. This helps to increase the device life.

If excessive candida overgrowth occurs treatment with antifungal medications may be considered as a preventive measure. Laboratory testing shows no negative influence on the function of the voice prosthesis when using the following antifungal medications: Nystatin, Fluconazole and Miconazol. Other drugs directly to or within close proximity of the voice prosthesis should be carefully assessed, since an unknown variety of chemical substances may influence the material.

2.4 Accessories

CAUTION: Only use genuine Provox accessories that are intended for use with Provox2. Other accessories may cause patient harm or product malfunction.

Accessories to facilitate insertion of the prosthesis (for clinician's use)

- **Provox Measure:** Indicates the length (prosthesis size) of the TE puncture.
- **Provox Dilator:** Is used to upsize smaller TE punctures to allow fitting of Provox voice prostheses or for blockage of a TE puncture. The dilator may also be used to temporarily block or stent the TE puncture.
- **Provox GuideWire:** For placement of a Provox2 or Provox Vega voice prosthesis in a freshly made TE puncture, or for retrograde replacement of a Provox indwelling voice prosthesis.

Accessories to maintain the device function (for patient's use)

- **Provox Brush / Provox Brush XL** are used to clean the lumen and valve of the prosthesis.
- **Provox Flush** allows for flushing the lumen of prosthesis with air or water.
- **Provox Plug** is a First Aid tool to temporarily block leakage through the prosthesis.

3. Adverse Events/Trouble Shooting Information

Prosthesis related

Problem: Leakage through the voice prosthesis. Leakage through the voice prosthesis is a normal development and may be an indication of the end of the device life.

Solution: Use a Provox Brush and/or Provox Flush to clean the voice prosthesis. If cleaning of the device does not resolve the leakage, replace the voice prosthesis.

Problem: Leakage through the voice prosthesis immediately after insertion.

Probable cause: Immediate leakage after voice prosthesis replacement can occur if the esophageal flange of the

prosthesis has not fully deployed into the lumen of the esophagus or if the prosthesis is too short. If the tip of the voice prosthesis inserter has not been placed far enough into the puncture during the insertion, the esophageal flange can unfold in the puncture tract instead of in the esophagus. Immediate leakage after voice prosthesis replacement can also occur if the voice prosthesis has been damaged during the insertion.

Solution: Rotate the voice prosthesis to check the fit. Perform an endoscopic examination to determine if the voice prosthesis is in the correct position and of the right length. Correct the position of the device. If the voice prosthesis is too short, replace with a longer device. If the leakage remains after repositioning, pull out the voice prosthesis with the inserter, or a non-toothed hemostat. Inspect the voice prosthesis. If the voice prosthesis has been damaged (e.g. the blue ring has been dislodged or the valve has been torn off), discard the prosthesis, otherwise reload and repeat the insertion procedure as described in section 2.1.

Handling related

Problem: It is difficult to push the prosthesis forward in the tube.

Probable cause: The esophageal flange has not folded completely forward into the loading tube.

Solution: Interrupt the insertion and pull the prosthesis back out of the loading tube. Inspect the voice prosthesis. If the voice prosthesis is damaged, discard the prosthesis! If the valve is undamaged reload and repeat the procedure as described in section 2.3.

Problem: The prosthesis has been completely inserted into the esophagus (accidental overshooting).

Solution: Leave the safety strap attached to the inserter and follow the instructions in section 2.4.

Puncture related

Aspiration of the prosthesis or other components – Immediate symptoms may include coughing, choking or wheezing.

Solution: Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If the patient can breathe, coughing may remove the foreign body. As with any other foreign body, complications from aspiration may cause obstruction or infection including pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary puncture and asthma.

Ingestion of the prosthesis or other components can cause symptoms that largely depend on size, location, degree of obstruction (if any).

Solution: The object usually passes spontaneously into the stomach and subsequently through the intestinal tract. If bowel obstruction, bleeding or perforation occurs, or the object fails to pass through the intestinal tract, surgical removal must be considered. If the device has remained in the esophagus it can be removed endoscopically.

Hemorrhage/Bleeding of the puncture – Slight bleeding from the edges of the TE-puncture may occur during replacement of the prosthesis and generally resolves spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage prior to placement or replacement of the prosthesis.

Infection and/or edema of the TE-puncture may increase the length of the puncture tract (e.g during radio therapy).

This may cause the prosthesis to be drawn inward and under the tracheal or esophageal mucosa. Inflammation or overgrowth of the esophageal mucosa may also cause the prosthesis to protrude from the puncture on the tracheal side. **Solution:** Temporarily replace the prosthesis with a longer one. Treatment with antibiotics with or without corticosteroids may also be considered. If the situation does not improve the prosthesis should be removed. Stenting the puncture with a catheter might be considered.

Granulation around the TE-puncture.

Solution: Electrical, chemical, or laser cauterization of the area of granulation may be considered.

Hypertrophic scarring around the puncture with bulging of the tracheal mucosa over the tracheal flange may occur if the prosthesis is relatively short.

Solution: The excess tissue may be removed by using a laser or a prosthesis with a longer shaft can be inserted.

Protrusion and subsequent extrusion of the prosthesis is sometimes observed during infection of the TE-puncture.

Solution: Removal of the prosthesis is required to avoid dislodgement into the trachea. Repuncture may be necessary if the puncture has closed.

Leakage around the prosthesis – Transient leakage may occur and can improve spontaneously.

Solution: The most common reason for persistent leakage around the voice prosthesis is that the device is too long. Replace the prosthesis with a shorter one. If the leakage is indicated by an enlarged puncture, place a Provox XtraFlange over the voice prosthesis' tracheal flange or replace it by a Provox Vega XtraSeal. Temporary removal of the prosthesis, a submucosal purse string suture or insertion of a cuffed tracheal cannula and/or nasogastric feeding tube may also be considered to reduce the size of the puncture. If the leakage persists, surgical closure of the puncture may be necessary.

Tracheal and/or esophageal tissue damage may occur if the prosthesis is too short, too long, or is pushed frequently against the esophageal wall by a tracheal cannula, stoma button, or the patient's finger. In patients undergoing radiotherapy this may happen more easily.

Solution: If the patient experiences soreness or pain in the area around the voice prosthesis, inspect the tissue around the TE puncture by endoscopy to avoid severe damage.

4. Date of Printing and User Assistance Information

The version number of this Instruction for Use and manufacturer contact information is located on the back cover.

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© Atos Medical AB, 2015.