

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	DD	Pontus Eklund - X-PONEKL	2020-04-21 - 14:00
Reviewed:	DD	Jon Berg - JONBER	2020-04-21 - 15:56
Approved:	DD	Fredrik Calais - FRECAL	2020-04-21 - 16:16
Released:	DD	Jon Berg - JONBER	2020-05-28 - 15:07

This document has been electronically signed by the persons above.

Product Information

Freevent® TubeBrush


Product description:

The Freevent TubeBrush is used for cleaning of tracheostomy tubes ex situ.

The Freevent TubeBrush is packed 10 pieces in a plastic bag.

It is available in five different models with outer diameter
 - 6 mm, 8 mm, 10 mm, 12 mm or 14 mm.

Document No: 10000031450 Edition: 02 Release date: 2020-05-28

Product Information

Document ID:	PF082-01-TechInfo	Edition:	02
Manufacturer:	Atos Medical AB Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden		
Classification: (EU) 2017/745	Class I, Rule 1		
Intended Use:	The Freevent TubeBrush is used for cleaning of tracheostomy tubes ex situ.		
CE Mark:	Yes, the devices are CE marked.		
GMDN code:	34883 (Airway device, cleaning brush, noninvasive).		
Sterilization:	Non-sterile.		
Raw material:	ABS, Stainless Steel, PBT and Cotton.		
Latex information:	Not manufactured with natural rubber latex		
Biological origin:	The device is not manufactured with materials derived from human or animal source.		
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.		
Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.		
Hazardous components:	None		
Expiration date:	3 years after manufacturing.		
Packaging:	10 brushes in each package except for Freevent TubeBrush set (1210) which includes 3 brushes (8, 10 and 12 mm) in each package.		

Devices under Basic UDI-DI: 7331791-GEN-A-000-0001-E9

REF	Name	UDI-DI
1205	Freevent TubeBrush Sz 6	7331791014833
1206	Freevent TubeBrush Sz 8	7331791014840
1207	Freevent TubeBrush Sz 10	7331791014857
1208	Freevent TubeBrush Sz 12	7331791014864
1209	Freevent TubeBrush Sz 14	7331791014871
1210	Freevent TubeBrush Set 1x8, 1x10, 1x12mm	7331791014888

Atos Medical AB Compatible products:

-